AIRSENSE A N A L Y T I C S



PEN3

Portable Electronic Nose

Intelligent chemical sensor for identification of gases and vapours

PEN3 is based on a metal-oxide gas sensor array. By using a specific dilution technique, the system is protected from overloading with substances. This extends the life time of the sensors and shortens the cycle times. This gasflow control makes it possible to stabilize the pattern under varying concentrations and analytical conditions.

Individual gases or gas mixtures are classified and recognized by their pattern. The unit then provides a quick and simple quantitative answer such as **good** or **bad**, **yes** or **no**

Qualitative answers can be developed using additional data and by comparing them with the database.

The variety of pattern recognition algorithms allows the system to be used in a broad range of applications - including laboraties and industries.

Process Control

- Process Control
- dosage of spices in food production
- supervision of industrial cleaning processes, fermentation processes
- dosage of artifi cial odour in natural gas
- production of polymer packaging material for food industry
- frying or roast process control



Laboratory equipment



Quality Control

Rancidity of oils, freshness of food, off odor in packaging materials, residual solvents in polymers, degradation of flavors, off odor in medicine, characterization of resins or aroma in beverages

Environmental & Safety Control

Odor in waste water purification plants or in compost plants (correlation with olfactometry), supervision of filters, solvents at workplace atmosphere, smouldering fires, identification of bacteria, leakage control or combustion emissions

Advantages

- Sensitive intelligent sensors
- Portable instrument with integrated display
- Small, fast, flexible and robust
- Online sampling technique
- Integrated dilution control
- Works with computer or in standalone mode
- Data logger and off-line analysis









PEN3

Portable Electronic Nose







Technical Data

Product Description

Sensor technology hot sensors, working temperature 350°C ... 500°C

Sensor array 10 different metal oxides single thick film sensors, optional adapted to application

Sensor response time typically: less than 1 second special fluidic connector

Sample flow adjustable, 10ml/min to 400ml/min

Flow system internal pumps, internal sample dilution system

Measurement time depending on the application from 4 seconds to some minutes

typically: 1 minute (20s measurements, 40s zero gas)

Sensitivity LOD 0.1 to 5ppm for gases and organic solvents, e.g. 0.1ppm H2S

Display graphical display
Dimensions 255 x 190 x 92 mm

Weight 2.1 kg

Environment Requirements

Temperature **typically:** 0°C to 45°C Humidity (relative) 5 % to 95 %, non-condensing

Power Requirements

Main power Power supply: 110 to 230VAC; max 30W or 12VDC (optional)

Communications

Computer interface USB port or serial RS-232 (optional)

System Requirements

Operating system Windows XP, Vista, Windows 7

Software WinMuster PEN for evaluation and analysis

Algorithm for analyse: Euclid, Correlation, Mahalanobis, DFA, PCA, LDA and PLS

Safety Class	Warranty
	Safety Class

Trap & Desorption Unit EDU Headspace-Autosampler	Compliant to EN292 Part1 & 2, EN294, EN61010-1. EN1050. EN60204-1. EN	12 months
neauspace-Autosampiei	55011 G1 CB. EN50270. EN61326	
	EN50270, EN61326	